Wednesday, September 12, 2018

Your Learning Goal:

After students learn the different parts of an electronic balance, they will practice using it to find the masses of 7 objects with 100% accuracy. Additionally, they will practice how to accurately predict the mass of common objects.

Table of Contents: Mass Mania – 4R & 4L

<u>Catalyst: (4L) Explain the image</u>

Homework: Observation & Inference DUE Tuesday/ Wednesday	<mark>Agenda:</mark> 1. Catalyst 2. Mass Mania 3. Reflection	
--	---	--

Table of Contents

<u>Date</u> 8/24/18 8/30/18 9/4/18 9/12/18 AssignmentPg #Marshmallow ChallengeObservation vs. InferenceRules of the RulerMass Mania *

1 L+ R 2 L+ R 3 L+ R 4 L + R

Catalyst:9/11/18Explain the image below



4L

<u>Mass Mania</u>

9/11/18

* What is the <u>electronic balance</u> used for?

Skip 6 lines

•Parts of the **Electronic Balance**:

- <u>Display Screen</u> *Skip 3 lines*
- <u>On/Off Button</u> <u>Skip 3 lines</u>
- Unit Button -



Catalyst: 4L



* <u>What is the electronic balance</u> <u>used for?</u>

- The <u>electronic balance</u> is used to <u>measure</u> the <u>mass</u> of an object. The <u>matter</u> that something contains.
- Weight is the pull of gravity on an object (and will change depending on location).
- The <u>unit</u> for <u>mass</u> is grams (g).



4R

* Parts of the <u>Electronic Balance</u>

<u>Unit Button</u> – Touch to change the unit to <u>grams</u> (g).



<u>Display Screen</u> – Shows you the amount of <u>mass</u>. Always wait until the numbers stop changing.

On/Off Button – Push once to turn on, hold down to turn off. Also makes the balance go back to zero.

Catalyst:

I think it <u>is/is not</u> possible to measure how much "stuff is inside an object without opening it because ...

LEAF:

<u>Mass Mania</u>

9/12/18

9/12/18

AR

* What is the <u>electronic balance</u> used for?

The <u>electronic balance</u> is used to measure the <u>mass</u> of an object.
The <u>unit</u> for <u>mass</u> is <u>grams</u> (g).

Parts of the <u>Electronic Balance</u>:
<u>Display Screen</u> - Shows you the amount of <u>mass</u>. Always wait until the numbers stop changing.

• <u>On/Off Button</u> – Push once to turn on, hold down to turn off. Also makes the <u>balance</u> go back to zero.

• <u>Unit Button</u> – check to make sure its in grams (g). Touching it will change the <u>unit</u>.

How to use the Electronic Balance



- Push the On/Off Button one time and the numbers should appear. The display should read 0.0 g.
- If the display has a different number like 0.1 or 0.2 hit the Zero button make the balance go back to zero.

What if my display doesn't read 0.0 g?

If your balance has 0.000 then your units of measurement are incorrect!



Press the Unit Button until the units change to grams (g). This should give you the correct units for measuring mass.

- Activity 1 Pre-lab:
 - Select four SMALL! single objects that you want to find the known <u>mass</u> for.
 - Write down the objects you selected before you place them on the balance. Remember, you are measuring <u>mass</u> so make sure the balance is measuring in the metric <u>units</u> of <u>grams (g)</u>.
 - Once the object is on the balance, *record* the <u>Mass</u> of the Single Object in <u>grams (g)</u>.

- Activity One
 - Select objects that you think will combine to equal the Target <u>Mass</u> required.
 - Write down the objects you selected before you place them on the balance. Remember, you are predicting so don't change your combination if it is wrong, just try to get more accurate as you go.
 - Once the objects are on the balance record the Measured Mass in grams (g).

- Activity Two
 - -For the next activity, *measure* the <u>mass</u> of three plastic Easter eggs and record the data in your table on the worksheet.
 - After you *measure* the <u>mass</u> of the eggs put the eggs in order from highest <u>mass</u> to least <u>mass</u>.
 - *Explain* which egg had the most <u>matter</u>
 based on what you've learned about <u>mass</u> so far.

- Activity Three
 - For the final activity you have a reading.
 Using what you have learned about <u>mass</u>, explain why there is a difference in the <u>mass</u> for the 3 eggs.

Mass vs. Weight





Lead: Where you state the topic of your paragraph. The mass of _____ and the mass of _____ were different because.....

Evidence: Observable and quantifiable data that a writer uses to support a claim. (Use your numerical data to support your lead)

Analysis/Warrant: Certain rules that connect evidence back to claims—how the evidence supports the claim. (What characteristics make the mass of these objects different?)

Finisher: Restating your claim in a new way to provide closure for your argument. (How are the mass and weight of your objects different? Why?)

Catalyst:

4

Measuring mass is different from measuring the weight of an object because...

<u>LEAF:</u> Describe how the mass and weight of two of the objects you measured are different and why.

The mass of _____ and the mass of _____ were different because.....

<u>Mass Mania</u>

9/11/18

9/11/18

* What is the <u>electronic balance</u> used for?

The <u>electronic balance</u> is used to measure the <u>mass</u> of an object.
The <u>unit</u> for <u>mass</u> is <u>grams</u> (g).

•Parts of the **Electronic Balance**:

•<u>Display Screen</u> - Shows you the amount of <u>mass</u>. Always wait until the numbers stop changing.

• <u>On/Off Button</u> – Push once to turn on, hold down to turn off. Also makes the <u>balance</u> go back to zero.

• <u>Unit Button</u> – DO NOT TOUCH, it will change the <u>unit</u>.